**Triple-Action Vacuum  
Relief Valve VRV-IIC**

**REF 4004103**

**1 UNIT**

**VACUUM RELIEF VALVE FOR  
LEFT VENTRICULAR DECOMPRESSION**

**DIRECTIONS:**

USE ASEPTIC TECHNIQUE.

- 1) Place in 1/4" line between left ventricle and suction pump with the "pump" end of the valve towards the pump. Arrows on the valve indicate flow direction, e.g., "Heart" toward the heart.

**PLEASE NOTE:** This valve will limit vacuum in the line when installed properly. If there is pressure in the line on the pump side, the internal check valve will block the flow of air towards the heart, and the air will be vented to the atmosphere through the annular opening of the side valve.


- 2) Do not obstruct valve opening.
- 3) Tighten all connections prior to use.

**NOTE:**

- 1) The valve may whistle if the dome is in a downward position
- 2) Cardioplegia line should be clamped off when vacuum line is activated. Vacuum line should be closed off when cardioplegia line is activated. Not clamping the cardioplegia line when activating vacuum line may cause positive pressure build up resulting in positive pressure relief to the atmosphere through the relief band.

For symbols glossary visit [www.questmedical.com/resources.aspx](http://www.questmedical.com/resources.aspx)

**STERILE EO** **Rx ONLY**

**STERILE** – in unopened or undamaged package.

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©2020 Quest Medical, Inc. All rights reserved.