
MiniGuard™ ARTERIAL SAFETY VALVE

REF 4008100/4008200 1 UNIT • APPAREIL • UNITÀ • EINHEIT • UNIDAD

The MiniGuard™ Arterial Safety Valve is intended to be used in an Extracorporeal Circuit during cardiopulmonary bypass procedures for the prevention of retrograde flow.

DIRECTIONS: USE ASEPTIC TECHNIQUE.

1. Remove from the package and discard the end caps.
2. Connect the valve to the 1/4" arterial line proximal to the arterial line filter (anywhere between the oxygenator outlet and the arterial filter) or anywhere retrograde flow could be harmful. The direction of flow is indicated on the valve by arrows. **WARNING:** The valve must be oriented properly with the direction of blood flow. Valve malfunction and possible damage to blood components will result if the valve is positioned incorrectly.
3. The valve must be primed in the forward flow direction.
4. When priming, hold the valve with the outlet end pointing up to facilitate the removal of air.
5. The valve and filter must be drained in the forward flow direction.

CAUTION: 1) This device is not intended for use at pressures above 10 psi nor for periods of use exceeding 6 hours. 2) If MiniGuard valve is placed proximal to arterial filter purge line, retrograde flow through the purge line may result. 3) Contents of un-opened, undamaged package are sterile/ non-pyrogenic. Do not resterilize or reuse.

FOLLOW INSTRUCTIONS CAREFULLY: The surgical team is solely responsible for determining the suitability of the valve relative to the CPB circuit, equipment, and surgical technique employed.

For symbols glossary visit www.questmedical.com/resources.aspx

QUEST Medical, Inc.

One Allentown Parkway, Allen, Texas 75002 USA
+1.972.390.9800 | +1.800.627.0226 | F +1.972.390.8465

STERILE EO

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Made in USA
904502 REV B
FEB 19

© 2019 Quest Medical, Inc. All rights reserved.

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