


REF 4170Y Approximate Priming Volume 0.09 mL



DESCRIPTION: Needle free accessory to a vascular access device

INTENDED USE:

Lyka[®] port needle free access device is intended for use as an accessory to a vascular access device (catheter) used in Hemodialysis or as an accessory to an I.V. Set for the administration or withdraw of fluids to a patient through a cannula or needle placed in the vein or artery.

CONTRAINDICATIONS: None known.

WARNINGS:

- 1) Intended for Single Patient Use only. **DO NOT REUSE.**
- 2) After use, this product may be a potential biohazard. Handle and discard in accordance with the accepted medical practice and applicable local, state, and federal laws and regulations.
- 3) For unattended medication delivery, only luer locking secondary sets and syringes should be used with the injection site.
- 4) Not for use longer than 7 days.

CAUTIONS:

- 1) Before using, examine the package carefully before opening to confirm its integrity and that the expiration date has not passed. The device is supplied in a sterile package and is non-pyrogenic. Do Not Use if package is damaged, opened or the expiration date has passed. Do not re-sterilize.
- 2) Access of the Lyka Port by syringe or I.V tubing should be straight on, not at an angle.
- 3) Do not use needles or blunt cannula with the Lyka Port.
- 4) Replace set per CDC guidelines or hospital protocol.
- 5) Flush injection site in accordance with hospital policy after or between injections.
- 6) To prevent retrograde flow after flushing injection site with a syringe, twist syringe counterclockwise and continue depressing syringe plunger until disengaged from site.

PROCEDURE - USE ASEPTIC TECHNIQUE:

- 1) Swab septum vigorously in circular motion with 70% alcohol antiseptic for 30 seconds and allow to air dry.
- 2) Prime in accordance with hospital policy
- 3) Attach male luer connector at distal end of syringe, catheter, or I.V. tubing into female luer of the Lyka Port by pushing in and twisting clockwise until resistance is felt. **Do not over tighten.** For slip luers, push syringe into the Lyka Port while twisting.
- 4) I.V. connection is ready for use. Deliver I.V. fluid or medication in usual manner.

DISCONNECTING I.V. LINES AND SYRINGES:

- 5) Disconnect by clamping line, then twisting male luer connector at the distal end of I.V. tubing or syringe counterclockwise until loose and apart. **DO NOT REMOVE THE LYKA PORT** from patient's catheter or tubing.
- 6) Flush the Lyka Port per hospital protocol.
- 7) To reconnect, repeat from step 1.
- 8) Replace per hospital protocol.

For symbols glossary visit www.questmedical.com/resources.aspx



STERILE EO

 **QUEST Medical, Inc.**

ONE ALLENTOWN PARKWAY
ALLEN, TEXAS 75002 USA
+1.972.390.9800 | +1.800.627.0226 | FAX +1.972.390.8465

Made in the U.S.A.
904583 REV A
JUL 19

© 2019 Quest Medical, Inc. All rights reserved.