



## 7 in. (18 cm) Extension Set

with 1 pressure activated checkvalve  
and gravity activated checkvalve

REF

**9250**

*Approximate priming volume 1.4 ml*

*1 UNIT*

### Intended Use:

For administration of intravenous fluids to a patient's vascular system utilizing needle-free components and I.V. manifold for multiple simultaneous intravenous therapy via gravity, syringe, or infusion pump.

### Procedure — Use Aseptic Technique:

1. Place pouch on surface and remove set.
2. Ensure all connections and caps are tightened as necessary.
3. Drip chamber sets: Close clamp, if applicable, on set. Prepare solution container. Insert spike and suspend container on I.V. pole. Squeeze and release drip chamber until it is approximately half full. Open clamp and allow fluid to fill tubing.
4. Non-drip chamber sets: Prime set by attaching main IV set to top most port on IV set and allow fluid to fill entire set.
5. Close clamp.
6. Connect set to patient's vascular access device and establish flow.

(OVER)

STERILE EO



Rx ONLY

LATEX FREE



For Warnings & Cautions  
Please consult:  
[www.questmedical.com](http://www.questmedical.com)

CE  
0 1 2 3

EC REP

M. Devices Group  
Marlborough House  
Riding Street  
Southport, PR8, IEW, UK  
Tel: +44 1704 544 944  
Fax: +44 1704 544 050



**QUEST** Medical, Inc.

An **Atrion** company [www.questmedical.com](http://www.questmedical.com)

One Allentown Parkway • Allen, TX 75002 • (800) 627-0226 • (972) 390-9800  
©2009 Quest Medical, Inc.

Made in U.S.A.  
904285 REV A  
SEP 09

#### **Secondary Infusions:**

Injection sites may have one-way injection ports in the following configurations:

- Pressure activated checkvalves (infusion via syringe or pump).
- Gravity check valves (infusion via gravity, syringe, or pump).

1. Prime a secondary set or syringe in the usual manner.
2. Remove cap.
3. Expel air as needed.
4. Attach secondary set and infuse desired amount of fluid.
5. Replace cap.

#### **IV Push:**

1. Prepare syringe.
2. Close clamp between injection site and vascular access device. Remove cap. Attach syringe. DO NOT allow air to enter line.
3. Open clamp downstream of injection site.
4. Push medication as prescribed and flush tubing.
5. Replace cap.