LaciCATH® Lacrimal Duct Balloon Catheter Instructions For Use

Description: The lacrimal duct catheter is a sterile, single-use, disposable balloon catheter consisting of a semi-flexible stainless steel core and nylon balloon tubing. The balloon is designed to inflate to a known diameter and length at the specific pressure. Markings are placed 10 mm and 15 mm proximal to the beginning of the working portion of the balloon, which indicate when the balloon is placed correctly in the lacrimal system. The overall length of the catheter is approximately 200 mm long. Pediatric LaciCATH lacrimal duct balloon catheters are available with 2mm or 3mm balloon diameter. For pediatric patients 30 months and older, the recommended catheter is the 3 mm model. For pediatric patients 12 months to 30 months of age, the recommended catheter is the 2mm model. The 3 mm balloon has a length of 15 mm and a folded profile of approximately 1.2 mm. The 2 mm balloon has a length of 13 mm and a folded profile of approximately 1.1 mm.

Intended Use Statement:

The LaciCATH Lacrimal Duct Balloon Catheter is a sterile, prescription device intended for single use during dilatation of an obstructed nasolacrimal duct system to correct epiphora.

Indications for Use Statement:

The LaciCATH Lacrimal Duct Balloon Catheter is a sterile, prescription device intended for single use during dilatation of an obstructed nasolacrimal duct system to correct epiphora as a result of the following:

2mm/3mm LaciCATH Pediatric catheter (intended for patients 12 months of age and older):

Congenital nasolacrimal duct obstruction
Dacryocystocele

Contraindications:

1. Lacrimal sac tumor
2. Dacryoliths

Warnings:

1. Intended for Single Patient Use only. DO NOT REUSE.
2. After use, this product may be a potential biohazard. Handle and discard in accordance with the accepted medical practice and applicable local, state, and federal laws and regulations.
3. Infection must be resolved with systemic antibiotics before surgery.
4. Deflate the balloon fully before withdrawing from the nasolacrimal duct.
5. The balloon must not be moved while inflated.

Cautions:

1. Do not exceed the recommended balloon inflation pressure.
2. Do not touch the balloon with instruments as it may puncture.
3. It is highly recommended that before using the lacrimal duct catheter in surgery for the first time, each surgeon reference a procedure overview sheet which can be downloaded from www.questmedical.com.

Possible Patient Complications:

The potential exists for serious complications including the following:

1. Excessive bleeding
2. Canaliculi trauma


Procedure:

1. Select size of catheter appropriate for patient.
2. Under general anesthesia, pass a Bowman probe down the nasolacrimal duct into the nose. Confirm the presence of the distal probe in the nose. Remove the probe.
3. CAUTION: Immediately before insertion of the catheter, remove the protective sleeve from the distal end. Lubricate with ophthalmic ointment if needed.
4. The balloon catheter is passed through the canaliculus and lacrimal sac into the nasolacrimal duct and pushed down the nasolacrimal duct to the nasal floor. Confirm the presence of the distal balloon catheter in the nose.
5. The balloon is gradually inflated to 8 atm (bars) for 90 seconds using a primed inflation device filled with 10 cc of sterile saline.
6. Deflate the balloon. If necessary, an additional inflation can be done by gradually inflating the balloon to 8 atm (bars) for 60 seconds and then deflate.
7. Pull the balloon proximally so that it lies within the lacrimal sac extending into the nasolacrimal duct. There are two markings on the proximal catheter shaft which are 10 mm and 15 mm from the beginning of the working portion of the balloon. The punctum is aligned halfway between the 10 mm mark and the balloon. The balloon can also be palpated delicately in the lacrimal location. The balloon is gradually inflated to 8 atm (bars) for 60 seconds.
8. Deflate the balloon. If necessary, an additional inflation can be done by gradually reinflating the balloon to 8 atm (bars) for 60 seconds.
9. Aspirate all residual fluid to fully deflate the balloon until a vacuum is drawn within the balloon as indicated on the inflator gauge.
10. Rotate the catheter clockwise to reduce the profile of the deflated balloon and withdraw the catheter.
11. Irrigate the lacrimal system with fluorescein and recover in the nose with a #5 or #8 flexible clear feeding tube used as a suction catheter to confirm patency.
12. If irrigation is not possible, another dilatation procedure can be performed using a new catheter.

LaciCATH® Lacrimal Duct Balloon Catheter

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For symbols glossary visit www.questmedical.com/resources.aspx

STERILE EO STERILE Rx ONLY CE 0 1 2 3

QUEST Medical, Inc.

One Allentown Parkway | Allen, Texas 75002 USA
+1.972.390.9800 | +1.800.627.0226 | F +1.972.390.8465

©2019 Quest Medical, Inc. LaciCATH is a registered trademark of Quest Medical, Inc.
The products displayed are protected by one or more patents worldwide, including US7169163, D679803, WO2004030720, ZL038227754 (China), and 259182 (Mexico).

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