## Lacrimal Duct Balloon Catheter

# LacriCATH® Instructions For Use

**Description:** The lacrimal duct catheter is a sterile, single-use, disposable balloon catheter consisting of a semi-flexible stainless steel core and nylon balloon tubing. The balloon is designed to inflate to a known diameter and length at the specific pressure. Markings are placed 10 mm and 15 mm proximal to the beginning of the working portion of the balloon, which indicate when the balloon is placed correctly in the lacrimal system. The overall length of the catheter is approximately 24 cm long with a working length of 40 mm. The catheter has a 5 mm balloon diameter with an approximate length of 8mm and an outer diameter deflated profile of approximately 1.6 mm.

**Intended Use:** The lacrimal duct catheter is a sterile, non-powered, prescription device intended for single use during dilation of an obstructed nasolacrimal duct system to correct epiphora.

**Indications For Use:** Adult LacriCATH catheters are indicated as a treatment for epiphora as a result of:

1. Functional or complete nasolacrimal duct obstruction
2. Dacryocystitis
3. Failed DCR (Dacryocystorhinostomy)

The 5 mm catheter is intended for use on adult patients.

### Contraindications:

1. Lacrimal sac tumor
2. Dacryoliths

### Warnings:

1. Intended for Single Patient Use only. DO NOT REUSE.
2. After use, this product may be a potential biohazard. Handle and discard in accordance with the accepted medical practice and applicable local, state, and federal laws and regulations.
3. Infection must be resolved with systemic antibiotics before surgery.
4. Deflate the balloon fully before withdrawing from the nasolacrimal duct.
5. The balloon must not be moved while inflated.

### Cautions:

1. Do not exceed the recommended balloon inflation pressure.
2. Do not touch the balloon with instruments as it may puncture.
3. It is highly recommended that before using the lacrimal duct catheter in surgery for the first time, each surgeon reference a procedure overview sheet which can be downloaded from [www.questmedical.com](http://www.questmedical.com).

### Possible Patient Complications:

The potential exists for serious complications including the following:

1. Excessive bleeding
2. Canaliculi trauma


### Procedure:

1. Dilate puncta.
2. Under local or general anesthesia, pass a #3 or #4 reinforced Bowman probe through the superior canaliculus into the lacrimal sac and push through the medial wall of the sac and the inferior part of the lacrimal fossa into the nose.
3. In the primary DCR, push the probe through three or four additional adjacent spots in the part of the inferior lacrimal fossa. Join these areas by moving the probe or reaching up the nose and inserting a nerve hook. In the failed DCR, the probe is brought through one spot of the obstructed ostium.
4. Remove the probe.
5. CAUTION: Immediately before insertion of the catheter, remove the protective sleeve from the distal end. Coat the distal end of the LacriCATH balloon catheter with ophthalmic ointment to facilitate insertion.
6. The balloon catheter is brought through the superior canaliculus, lacrimal sac, and out the opening in the medial lacrimal sac into the nose.
7. The balloon catheter is positioned with the aid of the exit markers on the distal working portion of the balloon, visible in the nose just past the opening in the lateral nasal wall. Markers are also placed at the proximal end of the balloon. The punctum should be aligned between the markings on the proximal catheter shaft, which are 10 mm and 15 mm from the beginning of the proximal working portion of the balloon.
8. Connect the catheter to a primed inflation device.
9. The balloon is gradually inflated to 8 atm (bars) for 90 seconds using the inflation device filled with 10 cc of sterile saline. Lateral traction can be placed on the catheter to prevent the balloon migration.
10. Deflate the balloon, then gradually reinflate to 8 atm (bars) for 60 seconds.
11. Visualize the ostium intranasally after dilation and confirm a good opening has been achieved. If not, the balloon should be repositioned and reinflated. CAUTION: Only move the balloon when deflated.
12. Aspirate all the residual fluid to fully deflate the balloon until a vacuum is drawn within the balloon as indicated on the inflator gauge.
13. Rotate the LacriCATH catheter clockwise to minimize the profile of the balloon and withdraw the catheter.
14. Be sure sac and ostium are dilated. A good opening should be visible after dilation.
15. Following the confirmation of a well-dilated opening in the ostium, place a STENTube® large diameter silicone tube through the opening. See STENTube instructions for use for complete placement instructions.
16. In the absence of contraindications, a course of systemic and oral corticosteroids can be used after dilation and later, after tube removal.

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For symbols glossary visit [www.questmedical.com/resources.aspx](http://www.questmedical.com/resources.aspx)

**STERILE** **EO**  **Rx ONLY** **CE**  
0 1 2 3

 **QUEST Medical, Inc.**

One Allentown Parkway | Allen, Texas 75002 USA  
+1.972.390.9800 | +1.800.627.0226 | F +1.972.390.8465

©2019 Quest Medical, Inc. LacriCATH is a registered trademark of Quest Medical, Inc.  
The products displayed are protected by one or more patents worldwide, including  
US7169163, D679803, WO2004030720, ZL038227754 (China), and 259182 (Mexico).

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