

Aortic Punches

CleanCut® & PerfectCut® Series

Create uniform aortic openings for vessel grafts

Produces uniform openings at the anastomotic site
Rotating cutting action requires less force
Self-aligning cutter mechanism
Six sizes and two lengths for a broad scope of surgical situations

Materials:

Medical-grade stainless steel
Thermoplastic construction

Available sterile and non-sterile from Quest Medical, Inc. and from all procedure kit suppliers

REF

SIZE

RCL28, RCC28	● 2.8 mm
RCL35, RCC35	● 3.5 mm
RCL40, RCC40, RCB40	● 4.0 mm
RCL45, RCC45, RCB45	● 4.5 mm
RCL48, RCC48, RCB48	● 4.8 mm
RCL50, RCC50, RCB50	● 5.0 mm
ALX33 (Lancet)	+ 3.3 mm

Punches 6 units/case | Lancet 25 units/case



custserv@questmedical.com

+800.627.0226 (option 1)
+1.972.390.9800 (v)
+1.972.390.7173 (f)

 **QUEST Medical, Inc.**
An **Atrion** company www.questmedical.com
© 2010 Quest Medical, Inc. CleanCut and PerfectCut are trademarks of Quest Medical, Inc. 2010-7