

OVER
15
YEARS
OF BETTER
OUTCOMES

LacriCATH[®] Procedure Kits







Your outcomes —
Our inspiration

Providing doctors the solution to chronic tearing

- Reduces trauma to the lacrimal system¹
- Eliminates placement of tubes
- High success rate¹
- Simultaneous bilateral inflation designed to reduce O.R. and anesthesia time

Sterile, single use procedure kits for the treatment of nasolacrimal duct obstruction

	Product Code	Description	Contents	Min Order Qty
2mm Dacryoplasty (DCP) 	DCP213-UNIT	Unilateral For patients under 30 months of age	(1) 2mm catheter (1) AQL [®] inflator	2
	DCP213-BIT	Bilateral For patients under 30 months of age	(2) 2mm catheters (1) AQL inflator (1) 4-way stopcock	1
3mm Dacryoplasty (DCP) 	DCP315-UNIT	Unilateral For patients over 30 months of age	(1) 3mm catheter (1) AQL inflator	2
	DCP315-BIT	Bilateral For patients over 30 months of age	(2) 3mm catheters (1) AQL inflator (1) 4-way stopcock	1
5mm Dacryocystorhinostomy (DCR) 	DCR508-UNIT	Unilateral For primary or revisional DCR	(1) 5mm catheter (1) STENTube [®] (1) AQL inflator	2
	DCR508-BIT	Bilateral For primary or revisional DCR	(2) 5mm catheters (2) STENTubes (1) AQL inflator	1
9mm Dacryocystorhinostomy (DCR) Developed by Bruce B. Becker, M.D. 	120DCR904-UNI	Unilateral 120° angle for easier lacrimal sac insertion	(1) 9mm catheter (1) STENTube (1) AQL inflator	2
	120DCR904-BI	Bilateral 120° angle for easier lacrimal sac insertion	(2) 9mm catheters (2) STENTubes (1) AQL inflator	1

www.questmedical.com

QUEST Medical, Inc.



© 2013 Quest Medical, Inc. LacriCATH and STENTube are registered trademarks of Quest Medical, Inc. QL is a registered trademark of Atrion Medical Products, Inc. The products displayed are protected by one or more patents worldwide, including US7169163, US6113567, D679803, WO2004030720, CN038227754 and PAA2005003261.

Ordering info

800.627.0226 • 972.390.9800

Fax: 972.390.7173

custserv@questmedical.com

Rx ONLY

2013-11

¹ Outcomes of 9mm Balloon-Assisted Endoscopic Dacryocystorhinostomy; Retrospective Review of 97 Cases. Silbert DI, Matta NS. ORBIT 2010;29(1):25-8

LacriCATH®

Product Components

Available Separately

Your outcomes —
Our inspiration

AQL1015 Atrion QL® Balloon Catheter Inflation Device



DCP/DCR accessories and instruments


Product Code	Description	Carton Qty
AQL1015	10cc inflation device for balloon catheter dilation. Lock design and gauge for accurate pressure control and monitoring. Measures pressures to 15 atm — sterile, single use	1
K12-01811	4-way stopcock for bilateral balloon DCP procedures, attaches to AQL inflation device for simultaneous bilateral inflation, included in 2mm and 3mm bilateral procedure kits — sterile, single use	10
LIS052	STENTube® large diameter silicone lacrimal intubation set for primary or revisional DCR — sterile, single use Developed by Bruce B. Becker, M.D.	2
LIS27T	Crawford-type silicone lacrimal intubation set, .025 in. tube diameter x 12.0 in. length, 4.5 in. probes olive tip ends — sterile, single use	3
Reusable CRH-K7	Crawford retrieving hook for use with STENTube and Crawford-type lacrimal intubation sets	1
Reusable RBP3-4	Reinforced Bowman probe for balloon DCR procedures	1
BLP1-2	Standard Bowman probe for balloon DCP procedures, size 1-2	1
BLP3-4	Standard Bowman probe for balloon DCP procedures, size 3-4	1
Reusable WLD024	Wilder dilator for balloon DCP and DCR procedures, 24mm	1
WLD019	Wilder dilator for balloon DCP and DCR procedures, 19mm	1
Reusable DNH001	Dandy nerve hook for balloon DCR procedures	1

K12-01811 Bilateral Stopcock



Lacrimal Intubation Sets

LIS052 STENTube®



LIS27T Crawford-type



CRH-K7 Crawford Hook



Bowman Probes

RBP3-4



BLP1-2



BLP3-4



Wilder Dilators

WLD024



WLD019



DNH001 Dandy Nerve Hook

