


Providing doctors the solution to chronic tearing:

- Reduces trauma to the lacrimal system
- Eliminates placement of tubes
- High success rate
- Simultaneous bilateral inflation reduces O.R. and anesthesia time





# LacriCATH<sup>®</sup> Procedure Kits

Sterile, single use procedure kits for the treatment of nasolacrimal duct obstruction



## LacriCATH Balloon Dacryoplasty (DCP) Procedure Kits • Unilateral and Bilateral

Kits include inflation device and catheter(s)

	PRODUCT CODE		CATHETERS/KIT	DESCRIPTION
	DCP213-UNI	Unilateral	1	2mm balloon diameter, 13mm length, balloon catheter for children less than 30 months of age.
	DCP213-BI	Bilateral	2	
	DCP315-UNI	Unilateral	1	3mm balloon diameter, 15mm length, balloon catheter for patients older than 30 months of age.
	DCP315-BI	Bilateral	2	

## LacriCATH Balloon Dacryocystorhinostomy (DCR) Procedure Kits • Unilateral and Bilateral

Kits include inflation device, catheter(s) and STENTube(s)<sup>®</sup>

	PRODUCT CODE		CATHETERS/KIT	DESCRIPTION
	DCR508-UNI	Unilateral	1	5mm balloon diameter, 8mm length, balloon catheter for primary DCR or revisional DCR.
	DCR508-BI	Bilateral	2	
	120DCR904-UNI	Unilateral	1	Transnasal balloon catheter for primary DCR or revisional DCR in adult patients. 9mm balloon diameter, 4mm balloon length, 120° angle for easier insertion in the lacrimal sac.
	120DCR904-BI	Bilateral	2	

# Product Components


Available Separately

QUEST Medical Ophthalmology Products



K12-01811

AQL1015



STENTube® LIS052

LIS27T

Crawford Hook CRH-K7

Bowman Probes RBP 3-4

BLP 1-2

BLP 3-4

Wilder Dilators WLD 024

WLD 019

Dandy Nerve Hook DNH001

## Atrion QL® Balloon Catheter Inflation Device

PRODUCT CODE	DESCRIPTION	QTY/KIT
AQL1015	One piece, plastic 10 cc inflation device for use with balloon dilation catheters. Lock design and gauge provide accurate control of inflation, monitoring balloon pressure and deflating balloon post-operatively. Measures pressures to 15 atm. Sterile, single-use.	1

## Bilateral Stopcock

PRODUCT CODE	DESCRIPTION	QTY/KIT
K12-01811	4-way stopcock used for bilateral balloon DCP procedures. Attaches to AQL1015 inflation device for simultaneous bilateral inflation. Standard in all 2mm and 3mm bilateral procedure kits. Sterile, single-use.	1

## Lacrimal Intubation Sets

PRODUCT CODE	DESCRIPTION	QTY/KIT
LIS052	<b>STENTube®</b> Lacrimal Intubation Set, developed by Bruce B. Becker, M.D. Large diameter silicone intubation tube.  Designed to provide ease of placement and enhanced stenting effect post DCR procedures. STENTube has two large diameter segments separated by a central thin segment to enhance placement and patient comfort.  Stainless steel probes are attached on tube ends and designed to resist separation. Sterile, single-use.	2
LIS27T	Silicone intubation set (Crawford-type). Silicone tube diameter: (.025 in) Tubing length: (12.0 in.) Probes: (4.5 in.) with olive tip ends. Sterile, single use.	3

## Instruments (Re-usable)

PRODUCT CODE	DESCRIPTION
CRH-K7	Crawford Retrieving Hook for use with STENTube® and Crawford-type lacrimal intubation sets.
RBP 3-4	Reinforced Bowman Probe for use with Balloon DCR procedures.
BLP 1-2	Bowman Probes for use with Balloon DCP procedures.
BLP 3-4	
WLD 024	Wilder Dilators for use with both Balloon DCP and DCR procedures.
WLD 019	
DNH001	Dandy Nerve Hook for use in Balloon DCR procedures.



**QUEST Medical, Inc.**  
An **Atrion** company

FOR ORDER INFORMATION: (800) 627-0226 • Fax: +1 972-390-7173  
One Allentown Parkway, Allen, TX 75002 • (972) 390-9800 • qmisales@questmedical.com • www.questmedical.com