


# Hands-free cardioplegia — Uninterrupted surgery

## Superior myocardial protection

With the Quest® MPS®2 Myocardial Protection System's Dual Delivery Line Set, there is no stopping to manage clamps or valves for any protection technique or strategy.



At your command —  
Antegrade, retrograde,  
or both.

**QUEST MPS<sup>®</sup>2**  
Myocardial Protection System  
*Protecting life with science*

The MPS2 provides independent control of potassium and cardioprotective drugs, providing the flexibility to deliver any cardioplegia protocol or strategy. The unique auto-pressure regulation capability maintains the optimal delivery pressure, maximizing flow distribution for the toughest ventricles, the hypertrophic heart. The Quest system allows surgeons and perfusionists to address the challenges of complex procedures, reduce blood utilization<sup>1</sup>, improve glucose management<sup>2</sup> and reduce costs.

**QUEST Medical, Inc.** [www.questmedical.com](http://www.questmedical.com) | 800.627.0226

1 Francesco Onorati, et al. "Polarizing" microplegia improves cardiac cycle efficiency after CABG for unstable angina.: Int J Cardiol; July 2012  
2 Bharat Datt MD, CCP, et al. Blood Only Cardioplegia with the Quest MPS. Canadian Perfusion Canadienne 2009; 19(1) April: 31-37  
©2014 Quest Medical, Inc. Quest and MPS are registered trademarks of Quest Medical, Inc. 2014-11

**R<sub>x</sub> ONLY**