



4007100 | 4007200

# RetroGuard<sup>®</sup>

## 3/8 in. Arterial Safety Valve

Optimum prevention of retrograde flow and related hazards

Automatically prevents retrograde arterial blood flow in the event of centrifugal pump failure


Large duckbill valve minimizes resistance to flow

Materials:

- Silicone duckbill valve
- Acetyl Butyl Styrene housing
- Latex free


Available sterile (4007100) and non-sterile (4007200) from Quest Medical, Inc. and from all procedure kit suppliers

RetroGuard<sup>®</sup>



# RetroGuard<sup>®</sup>

*Automatically prevents retrograde arterial blood flow in the event of centrifugal pump failure*



custserv@questmedical.com

+800.627.0226 (option 1)  
+1.972.390.9800 (v)  
+1.972.390.7173 (f)

 **QUEST Medical, Inc.**  
An **Atrion** company [www.questmedical.com](http://www.questmedical.com)  
© 2010 Quest Medical, Inc. RetroGuard is a trademark of Quest Medical, Inc. 2009-14