



4007100 | 4007200

# RetroGuard<sup>®</sup>

## 3/8 in. Arterial Safety Valve

Optimum prevention of retrograde flow and related hazards


Automatically prevents retrograde arterial blood flow in the event of centrifugal pump failure

Large duckbill valve minimizes resistance to flow

Materials:


- Silicone duckbill valve
- Acetyl Butyl Styrene housing
- Latex free

Available sterile (4007100) and non-sterile (4007200) from Quest Medical, Inc. and from all procedure kit suppliers



# RetroGuard<sup>®</sup>

*Automatically prevents retrograde arterial blood flow in the event of centrifugal pump failure*



 **QUEST Medical, Inc.** [www.questmedical.com](http://www.questmedical.com)

© 2014 Quest Medical, Inc. RetroGuard is a trademark of Quest Medical, Inc. 2009-14

[custserv@questmedical.com](mailto:custserv@questmedical.com)

**Rx ONLY**

+800.627.0226 (option 1) | +972.390.9800 (v) | +972.390.7173 (f)