

# INTRODUCING

## The All-New **MiniGuard™**

1/4 IN. ARTERIAL  
SAFETY VALVE

The new 1/4 in.  
MiniGuard™ with the  
3/8 in. RetroGuard®

Optimum prevention  
of retrograde flow

The MiniGuard™ is intended to be used in an Extracorporeal Circuit during cardiopulmonary bypass procedures for the prevention of retrograde flow and can be incorporated in several different applications where 1/4 inch tubing is needed.


Priming Volume: 6 mL

Opening Pressure:  $\leq .5$  psi

Flow Rate: Between 35 to 3500 mL



Automatically prevents  
retrograde arterial blood flow



 **QUEST Medical, Inc.**

custserv@questmedical.com | +1.800.627.0226  
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com

Rx ONLY



©2019 Quest Medical, Inc. RetroGuard is a registered trademark of Quest Medical, Inc. 2018-08