



MULTIPORT® MANIFOLD


SIMPLIFY & PROTECT (9520)

Six Ports, One Flow

Give your patients the quality care they deserve with Q2® split septum, needleless Multiport® Manifold. In-line gravity check-valves provide one-way flow without manipulation of stopcock handles. Quest's Q2 fluid delivery systems offer convenient designs, configured to simplify the delivery of medications.

MultiPort® 9520 is designed to mitigate risk and increase the efficiency of patient care.

- Closed system designed to reduce the chance of contamination, maintaining a sterile environment.¹
- The Q2 fluid delivery systems are configured to simplify the delivery of medications.
- Straight, direct fluid pathways for optimum flushing.
- Low residual volume/dead space.
- Neutral fluid displacement when attaching and detaching.
- Q2 check-valves ensure one-way fluid delivery of medication to the patient, preventing retrograde flow.
- Q2 MultiPort Manifold Extension Set offers 6 ports, 7 gravity-activated check valves and on/off pinch clamp



Exclusively Distributed by:



997 Horan Drive | Fenton, MO 63026
314.961.5786 | 800.969.6331 | cs@progressivemedinc.com
www.progressivemedinc.com



custserv@questmedical.com | +1.800.627.0226
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com

Rx ONLY



¹ "Closed catheter access systems are associated with fewer CRBSIs than open systems and should be used preferentially." - The Centers for Disease Control and Prevention of Intravascular Catheter-Related Infections, pg. 55