


# LacriCATH<sup>®</sup>

## Crawford Lacrimal Intubation Set (LIS27T)

Crawford Tubes are a highly effective tool designed to treat canalicular system stenosis or nasolacrimal duct obstruction utilizing non-traumatic stainless steel probes with olive tips.<sup>1</sup>

Type of tubing	Medical-Grade Silicone
Diameter	.025 in. (0.6 mm)
Length	12 in. (150 mm)
Sterile, disposable	
Box of 3	

For competitive pricing or to request a sample, please contact your LacriCATH sales representative today.



 **QUEST** Medical, Inc.

1 Tong NX, Zhao YY, Jin XM. Use of the Crawford tube for symptomatic epiphora without nasolacrimal obstruction. Int J Phththalmol 2016;9(2):282-285

© 2019 Quest Medical, Inc. 2019-23



[custserv@questmedical.com](mailto:custserv@questmedical.com) | +1.800.627.0226 | [questmedical.com](http://questmedical.com)