

INTRODUCING

EXOFLEX™



The manifold that morphs with your clinical demands

ExoFlex™ with Q2® Technology is a new manifold designed to span departments and accommodate the need of being **rigid or flexible - in seconds.**

- Adaptable between departments and clinicians
- Reduce risk of HA-BSI with the Q2 swabable, split septum valve¹
- Low residual volume/dead space
- Non-hemolytic Q2 swabable valve for repeated access²
- Integrated check valves ensure one-way infusion to the patient; preventing retrograde flow
- High flow rates for bolus needs
- Compatible with ISO 594 / 80369-7 disinfectant caps containing 70% IPA²
- ExoFlex decreases line manipulation reducing associated errors and infection risks



ExoFlex. Spanning departments without breaking the line.

Contact us for a sample or for more details.

Exclusively Distributed by:



997 Horan Drive | Fenton, MO 63026 | 314.961.5786 | 800.969.6331
cs@progressivemedinc.com | www.progressivemedinc.com



R_x ONLY

QUEST Medical, Inc.

custserv@questmedical.com | +1.800.627.0226
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com

¹ CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services; 2 Data on file
© 2019 Quest Medical, Inc. 2019-31C