

Infusion Therapy Products



Distributed by:

Product Description Key

Product listings contain one or more of the following descriptive icons in the product description key below. Some icons represent ideal usage but does not limit any product to an implied use.



Q2® Split Septum Technology



Low Volume (<3 mL)



Aspiration Port



CE marked



Approximate Priming Volume



Lead Time:
Thirty working days



Quest Medical, Inc. located in Allen, Texas engineers and manufactures a full line of infusion therapy products, including rigid and flexible manifolds, removable spine components, and intravenous (IV) administration sets in high or low volume configurations. Quest's R&D department continues to develop and improve offerings for the Infusion Therapy product lines.

Quest has revolutionized split septum technology with Q2® needleless technology - enhanced safety features that assist in reducing infection rates. Integrated into most Quest products, the Q2 straight unobstructed fluid path provides high flow rates and is engineered to suppress thrombotic occlusion risks caused by disconnecting syringes or secondary sets. Q2 needleless connectors enable multiple activations and maintain a sterile, closed system for increased patient safety and improved outcomes.

Our product lines include:
the New ExoFlex™ manifolds, Q2 luer access ports, extension sets,
and IV administration sets, and MultiPort® manifolds.

Quest products serve and support several areas of practice
including: Anesthesia, OR, ICU, Oncology, CCU, Cardiology,
Stem Cell, and other specialty departments.



One Allentown Parkway | Allen, Texas 75002-4211 USA
+1.800.627.0226 | +1.972.390.9800 | Fax: 972.390.7173
custserv@questmedical.com | questmedical.com

The products listed within this catalog are covered by one or more pending or issued patents including:
US D569391 S1, US D557281 S1, US D557282 S1, US D557283 S1, US D557285 S1, US D570880 S1, US D570373 S1

CONTENTS

NEEDLELESS CONNECTORS

3150	3
3160	3

PRECISION DELIVERY®

PD60	4
PD60F	4
PDM60	5
PDM60F	5
PDMV60	5
PDMV60F	5
PDV60	5
PDV60F	5
PDV06	5
PDV06F	5
PDT05	6
PDT05F	6
PDT05DF	6
PDTV05	6
PDTV05F	6
PDTV05DF	6

EXOFLEX™

95786	7
95787	7
95770	7
95789	7
95781	8
95782	8
95788	8
95779	8
95778	8

MANIFOLDS

95712	9
95711	9
9520	9
95710	9
95704	9
95702	9
95713	10
95707	10

EXTENSION SETS

9525A	10
9219	10
95701	10
95725S	10
95902	11
95907	11
95906	11
95714	11

IV ADMINISTRATION SETS

9527B	12
95780	12
22-201-V	12
9526B	12



Quest Split Septum Needless Connectors

“ When needless systems are used, a split septum valve may be preferred over some mechanical valves due to increased risk of infection with the mechanical valves.”

– CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services



Septum surface easily disinfected
Smooth, gapless surface

Septum seal reduces risk of biofilm formation
Housing tight between septum

Fluid pathway facilitates optimum flushing
Straight and direct

Heparin-free flushing solution
Uses saline to reduce the risk of heparin-induced thrombocytopenia


No moving parts reduces risk of HA-BSI contamination
No internal mechanism

Needless Connectors

Quest Medical, Inc. has revolutionized split septum technology with enhanced safety features that assist in reducing infection rates - the **Q2® Technology**. Integrated into most Quest products, the Q2 needless connector's straight, unobstructed fluid path provides high flow rates and is engineered to suppress thrombotic occlusion risks caused by disconnecting syringes and other sets. Q2 needless connectors enable multiple activations and maintain a sterile, closed system for increased patient safety and improved outcomes.



3150




Q2® swabable, needless connector

TUBING ID	APPROX. LENGTH	APV	CASE QTY
N/A	N/A	0.09 mL	200



3160



Q2® swabable, needless vial adapter

TUBING ID	APPROX. LENGTH	APV	CASE QTY
N/A	N/A	0.06 mL	50






Precision Delivery Extension Sets

Quest Medical's Precision Delivery® Extension Sets are intentionally engineered for safe and precise infusions. Specialized componentry minimizes dilution and priming volume while delivering all solutions to neonatal and pediatric patients. Our offerings include ultra-microbore tubing, pediatric-specific filters, Q2® Split-Septum Needleless Valve Technology, and a variety of configuration options for ultimate control of infusions.

PD60




60 in. ext with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	61 in. (155cm)	1.4 mL	100



PD60F




60 in. ext., 1.2 micron filter with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	62 in. (158cm)	1.7 mL	100



PDM60




60 in. ext., ultra microbore tubing, with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	61 in. (155cm)	0.5 mL	100



PDM60F




60 in. ext., ultra microbore tubing, 0.2 micron filter, with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	62 in. (158cm)	0.6 mL	100



PDMV60




60 in. ext., needleless connector, ultra microbore tubing with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	61 in. (155cm)	0.6 mL	100



PDMV60F



60 in. ext., needleless connector, ultra microbore tubing, 0.2 micron filter with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	61 in. (155cm)	0.8 mL	100



PDV60




60 in. ext, needleless connector with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	61 in. (155cm)	1.6 mL	100



PDV60F



60 in ext., needleless connector with, 1.2 micron filter with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	62 in. (158cm)	1.7 mL	100



PDV06




6 in. ext., needleless connector, with sliding clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.75 in. (20cm)	0.4 mL	50



PDV06F



6in. ext, needleless connector, ultra microbore tubing, 1.2 micron filter with slide clamp


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	8.5 in. (22cm)	0.7 mL	50



PRECISION DELIVERY EXTENSION SETS

continued

PDT05




5 in. ext., yellow slide clamp,
blue slide clamp, white slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7 in. (18cm)	0.2 mL (W)	50
		0.2 mL (B)	
		0.2 mL (Y)	



PDT05F




60 in. ext., ultra microbore tubing,
0.2 micron filter, with slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.25 in. (22cm)	0.2 mL (W)	50
		0.2 mL (B)	
		0.4 mL (Y)	



PDT05DF




5 in. ext., 0.2 micron filter, 1.2 micron
filter, yellow slide clamp, blue slide
clamp, white slide clamp and xv

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.25 in. (18cm)	0.4 mL (W)	50
		0.2 mL (B)	
		0.4 mL (Y)	



PDTV05




5 in. ext., needleless connector (3),
yellow slide clamp, blue slide clamp,
white slide clamp and 3-in-1 male luer

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.25 in. (18cm)	0.25 mL (W)	50
		0.25 mL (B)	
		0.25 mL (Y)	



PDTV05F




5 in. ext., needleless connector (3), 0.2
Micron Filter, yellow slide clamp, blue
slide clamp, white slide clamp and
3-in-1 male luer

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.25 in. (18cm)	0.25 mL (W)	50
		0.25 mL (B)	
		0.50 mL (Y)	



PDTV05DF



Needless Connector (3), 0.2 micron
filter, 1.2 Micron Filter, yellow slide
clamp, blue slide clamp, white
slide clamp and 3-in-1 male luer

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.25 in. (18cm)	0.50 mL (W)	50
		0.25 mL (B)	
		0.50 mL (Y)	



EXOFLEX™




The manifold that morphs with your clinical demands

Quest Medical has expanded the flexible manifold line with the new ExoFlex™ product line with a removable spine. Infection reduction was an inspiration in the development our new ExoFlex line; intended to address the needs of multiple departments while preventing a break in the line.

Features & Benefits of ExoFlex™

- Adaptable between departments and clinicians with a quick and easy to remove spine
- Reduce risk of HA-BSI with the Q2® split septum valve¹
- Integrated check valves ensure one-way infusion to the patient; preventing retrograde flow¹
- High flow rates for bolus delivery needs
- Reduce soft costs as the product is assembled, sterilized and ready for use

95786




ExoFlex™

Needleless connector with integrated check valve, (4) double T-site needleless connectors with integrated check valves and removable spines

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	11.25 in. (29cm)	1.75 mL	50



95787




ExoFlex™

Needleless connector with integrated check valve, (6) T-site needleless connectors with integrated check valves and removable spines

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	13.75 in. (37.5cm)	1.65 mL	50



95770




ExoFlex™

Needleless connector with integrated check valve, (2) double T-site needleless connectors with integrated check valves and removable spine

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	9 in. (22cm)	1.2 mL	50



95789




ExoFlex™

Needleless connector with integrated check valve, (2) double T-site needleless connectors with integrated check valves and removable spine, (1) Y-site connector

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	53 in. (135cm)	6.6 mL	50



95781




ExoFlex™

(3) T-site needleless connectors with integrated check valves and removable spine

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	10.5 in. (26.75cm)	1.75 mL	50



95782



ExoFlex™

Needleless connector with integrated check valve, (3) T-site needleless connectors with integrated check valves and removable spine

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	10.5 in. (26.75cm)	1.17 mL	50



95788



ExoFlex™

Needleless connector with integrated check valve, (3) T-site needleless connectors with integrated check valves and removable spine, (1) Y-site connector

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	55 in. (140cm)	7.3 mL	50



95779




ExoFlex™

Needleless connector with integrated check valve, (2) T-site needleless connectors with integrated check valves and removable spine, (1) Y-site connector

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	53 in. (135cm)	6.7 mL	50



95778



ExoFlex™

Needleless connector with integrated check valve, (2) T-site needleless connectors with integrated check valves and removable spine


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	9 in. (22cm)	1.2 mL	50



Manifolds

Quest Medical offers a variety of manifold options to assist with safe and simple solutions for fluid management in both rigid and flexible configurations. Incorporated into most of our infusion therapy products is the Q2 technology - a split septum, needleless connector. Q2 connectors offer enhanced safety features that assist in reducing infection rates, have a straight unobstructed fluid path that provide high flow rates and allow for multiple activations while maintaining a sterile, closed system for increased patient safety and improved outcomes.

95712




Q2® MultiPort®

(6) needleless connectors with integrated check valves,
(2) additional T-site needleless connectors

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	10 in. (25cm)	2.2 mL	50



95711




Q2® Double T-Site

Needleless connector with integrated check valve, (3) double T-site needleless connectors with integrated check valves, (2) slide clamps and additional T-site needleless connector

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	8.5 in. (22cm)	1.5 mL	50



9520




Q2® MultiPort®

(6) needleless connectors with integrated check valves

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	10 in. (25cm)	2.0 mL	50



95710




Q2® T-Site

Needleless connector with integrated check valve, (6) T-site needleless connectors with integrated check valves

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	14.5 in. (37cm)	1.0 mL	50



95704




Q2® Double T-Site

Needleless connector with integrated check valve, (2) double T-site needleless connectors with integrated check valves

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	7 in. (18cm)	1.1 mL	50



95702



Q2® T-Site


Needleless connector with integrated check valve, (3) T-site needleless connectors with integrated check valves

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	10 in. (25cm)	0.9 mL	50



MANIFOLDS *continued*

95713




Q2® T-Site

- (1) needleless connector with integrated check valve,
- (3) T-site needleless connectors

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	10 in. (25 cm)	0.8 mL	50



95707



Q2® T-Site

- (1) needleless connector with integrated check valve,
- (2) T-site needleless connectors


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	4 in. (10 cm)	0.6 mL	50



Extension sets

Quest's extension sets focus on safe and accurate infusion therapy delivery with minimal priming volumes and high flow rates. Our offerings include various tubing sizes and component options to ensure your patients are kept safe.

9525A



Q2® MultiPort®

- (6) needleless connectors with integrated check valves, high-flow 0.22 micron filter,
- (3) pinch clamps, in-line check valve and (3) Y-sites

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	100 in. (254cm)	15.5 mL	48



9219



Female Y-Extension Set

- (1) female luer pressure-activated check valve,
- (1) needleless connector with integrated check valve

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.105 in.	42 in. (107cm)	6.2 mL	50



95701




Q2® T-Site

- Needleless connector with integrated check valve, (1) T-site needleless connector with integrated check valve

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.06 in.	7 in. (18cm)	0.66 mL	50



95725S



Multilaminate IV Extension Set

- with 24 in. perforated sleeve

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.045 in.	25 ft (7.6m)	10 mL	50



EXTENSION SETS *continued*

95902




High flow filter Extension Set
High-flow 0.22 micron filter, pinch clamp and Y-site

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	17 in. (43.2cm)	4.5 mL	100



**LOW PRESSURE, POWER
INJECTABLE SETS** *up to 325 psi*

95907



Q2® Low Pressure Power Injection Extension Set
Needleless T-connector extension with needleless side port

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.045 in.	6 in. (15cm)	0.4 mL	100



95906



Q2® Low Pressure Power Injection Extension Set
Needleless T-connector extension

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.045 in.	6 in. (15cm)	0.3 mL	100



95714



Q2® Low Pressure Power Injection Extension Set
with detachable needleless connector


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.10 in.	7 in. (18cm)	1.0 mL	100



IV ADMINISTRATION Sets

Quest’s IV administration sets focus on safe and accurate infusion therapy by incorporating the use of split-septum neutral valves and preventing retrograde flow. Our variety of offerings are engineered to meet the demands of various departments in your facility including Anesthesia, Oncology and Critical Care.

9527B




Q2® MultiPort®
10-drop/mL admin set with (6) needleless connectors with integrated check valves, roller clamp, high-flow 0.22 micron filter, check valve and (3) Y-sites

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	115 in. (292cm)	19.0 mL	48



95780




ExoFlex™
10-drop/mL admin set with check valve, Y-site, roller clamp, (3) T-site needleless connectors with integrated check valves with removable spine and removable extension set with (2) Y-sites

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	117 in. (297.18cm)	22 mL	50




22-201-V



20-drop/mL vented secondary set with ball shut-off valve and roller clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.01 in.	75 in. (190cm)	15.5 mL	48

9526B



Q2® MultiPort®
60-drop/mL admin set with (6) needleless connectors with integrated check valves, roller clamp, high-flow 0.22 micron filter, check valve and (3) Y-sites

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.105 in.	115 in. (292cm)	19.0 mL	48





custserv@questmedical.com | +1.800.627.0226
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com

© 2020 Quest Medical, Inc. Q2, Checkmate, and Multiport are registered trademarks of Quest Medical, Inc. 2020-21
1. CardioMed Supplies, Inc. distributes in Canada (all provinces) and Bermuda.



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

