



Infusion Therapy Products



QUEST Medical, Inc.



Quest Medical, Inc. is committed to partnering with clinicians worldwide to produce better patient outcomes. We engineer infusion therapy products for superior infection control, organization, versatility and ease of use. Our innovative and high quality IV delivery sets are used to supply vital fluids and medications to patients in every department of your hospital or clinic.

Quest Medical aims to create efficiencies with each partnership through:

- Inventory Management**
- Reliable Client Services**
- Next Day Shipping**

Our representatives are subject matter experts who will assist you in standardizing products across departments, resulting in cost savings, reduced waste and faster set-up time.

To locate a sales specialist in your area, contact **800-672-0226** or email **custserv@questmedical.com**.

Icon Description Key



Contains Q2® Valve Technology



Low Volume (<3mL)




Contains Aspiration Port



European Conformity



Lead Time: 30 working days



BLOOD ADMINISTRATION SETS

Y115
YV115
YV130

5

IV ADMINISTRATION SETS

ASV10125
ASV10115
ASV10130
ASV60125

6

ASV101Y
ASV102Y
ASV103Y
95780
22-201-V

7

MICROBORE ADMINISTRATION SETS

ASV170F
ASV370
ASV470

8

MANIFOLDS

95786
95787
95770
95789

9

95781
95782
95791
95788
95779
95778
95712

10

95711
9520
95710
95704
95702
95707

11

HIGH-FLOW STOPCOCKS

ST01
ST02
ST03
STV01
STV02
STV03

12

HFSTV01
HFSTV02
HFSTV03

13

PRECISION DELIVERY® EXTENSION SETS

PD06F
PDM06F
PDMV06F
PDV06
PDV06F

14

PD60
PD60F
PDM60
PDM60F
PDMV60
PDMV60F
PDV60

15

PDV60F
PDBV05
PDB05F
PDB05DF
PDBV05F
PDBV05DF

16

PDT05
PDT05F
PDT05DF
PDTV05
PDTV05F
PDTV05DF

17

LIGHT-PROTECTIVE, AMBER EXTENSION SETS

PDV06A
PDV60A
PDMV60FA
PDV60FA
PDTV05A
PDTV05FA

18

PDTV05DFA

19

EXTENSION SETS

95701
9219
95902
95725S
95724

20

95725
95726
95727
HF0135

21

LOW-PRESSURE, POWER INJECTION EXTENSION SETS

95714
95906
95907

Q2® NEEDLELESS CONNECTORS

3150
3150HF
3160

22

Our Evolving Product Lines

Quest Medical, Inc. was founded on the premise of developing medical devices that would have a positive impact on patient care. Since 1979, we have been dedicated to providing the highest quality products to our three focused markets: Infusion Therapy, Cardiovascular and Ophthalmology.


Quest's evolving Infusion Therapy Line features rigid and flexible manifolds with removable components and intravenous (IV) administration sets with high or low volume configurations to fit all patients' needs. Our products incorporate the following trademarks: Q2® Valve Technology (Needleless Connectors/NCs), Precision Delivery® Extension Sets, ExoFlex® Manifolds, MultiPort® Manifolds and Qlinch™ Securement Bands.



Integrated into most of Quest's Infusion Products is our Q2® Valve Technology. Gaining access to each patient's bloodstream via Infusion Therapy can increase risk of infection. Minimize that risk by incorporating a simple, closed-system design into your line setup. Quest's Q2® Valve Technology is void of internal componentry, reducing the chance of contamination.


“When needleless systems are used, a split septum valve may be preferred over some mechanical valves due to increased risk of infection with the mechanical valves.”

- CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services




Blood Administration Sets


Y115	TUBING ID	APPROX. LENGTH	APV	CASE QTY
(2) high-flow vented spikes, (3) roller clamps, 200 micron filter chamber, large volume cylinder hand pump, slide clamp, (3) Y-sites + NCs, 32 in. detachable high flow 4-way 2-gang stopcock extension with pinch clamp, and detachable 6 in. extension with pinch clamp	.12 in.	119 in. (302 cm)	76 mL	50



YV115	TUBING ID	APPROX. LENGTH	APV	CASE QTY
(2) high-flow vented spikes, (3) roller clamps, 200 micron filter chamber, large volume cylinder hand pump, slide clamp and 35 in. detachable high flow 4-way 2-gang stopcock + Q2® NCs extension	.12 in.	114 in. (290 cm)	74 mL	50



YV130	TUBING ID	APPROX. LENGTH	APV	CASE QTY
(2) high-flow vented spikes, (3) roller clamps, 200 micron filter chamber, large volume cylinder hand pump, (2) slide clamps, (2) Y-sites + Q2® NCs, 45 in. detachable high flow 4-way 3-gang stopcock + Q2® NCs extension, and detachable 10 in. extension	.12 in.	129 in. (329 cm)	75 mL	50



IV ADMINISTRATION SETS



Quest's IV Administration Sets assist with accurate delivery of medication. Accessing the patient's bloodstream when administering drugs may increase the risk of infection. With the incorporation of Q2® Needleless Connectors, Quest's IV Administration Sets are engineered with a closed system to promote patient safety. Multiple tubing ID options and detachable extensions provide a variety of configurations, which meet the demands of various departments in your facility including Anesthesia, Oncology and Critical Care.



ASV10125


Vented 10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, 40 in. detachable high-flow 4-way 2-gang stopcock extension + Q2® NCs and pinch clamp

TUBING ID
0.10 in.

APPROX. LENGTH
129 in. (327cm)

APV
28.6 mL

CASE QTY
50



ASV10115


Vented 10-drop/mL administration set with check valve, Y-site with Q2® NC, roller clamp, 22 in. detachable high-flow 4-way stopcock extension+ Q2® NCs and detachable 10 in. high-flow 4-way 2-gang stopcock extension + Q2® NCs

TUBING ID
0.10 in.

APPROX. LENGTH
115 in. (292cm)

APV
27.9 mL

CASE QTY
50



ASV10130


Vented 10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, detachable 46 in. high-flow 4-way 3-gang stopcock extension + Q2® NCs and detachable 6 in. extension with Q2® NC and pinch clamp

TUBING ID
0.10 in.

APPROX. LENGTH
132 in. (335cm)

APV
30.4 mL

CASE QTY
50



ASV60125


Vented 60-drop/mL administration set with check valve, (4) Y-sites with Q2® NCs, (2) pinch clamps, roller clamp, 35 in. detachable high-flow 4-way stopcock extension + Q2® NC and detachable 6.5 in. extension with Q2® NC

TUBING ID
0.10 in.

APPROX. LENGTH
126 in. (320cm)

APV
29.1 mL


CASE QTY
50



IV ADMINISTRATION SETS


ASV101Y

10-drop/mL administration set with check valve, (1) Y-site with Q2® NC, roller clamp, (1) T-site with Q2® NC + integrated check valve and pinch clamp




ASV102Y

10-drop/mL administration set with check valve, (2) Y-sites with Q2® NCs, roller clamp, (1) T-site with Q2® NC + integrated check valve and pinch clamp




ASV103Y

10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, (2) T-sites with Q2® NCs + integrated check valve and pinch clamp



95780

10-drop/mL administration set with check valve, (3) Y-sites + Q2® NC, roller clamp, ExoFlex™ Manifold with (3) T-sites with Q2® NCs + integrated check valves + removable spine and detachable 45 in. extension set with pinch clamp



22-201-V

20-drop/mL vented secondary set with ball shut-off valve and roller clamp




IV ADMINISTRATION SETS




Microbore Administration Sets


ASV170F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore administration set with high-flow spike, pressure activated valve, pinch clamp, Q2® NC syringe port, check valve, slide clamp and 0.2 micron filter	0.10/0.04 in.	107 in. (272cm)	2.8 mL*	50



ASV370	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore trifurcated administration set with high-flow spike, pressure activated valve, pinch clamp, (3) Q2® NC syringe ports, check valve and (3) slide clamps	0.10/0.04 in.	103 in. (261cm)	0.3 mL (W) 0.3 mL (Y) 2.8 mL (B)*	50



ASV470	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore quadfurcated administration set with high-flow spike, pressure activated valve, pinch clamp, (4) Q2® NC syringe ports, check valve and (4) slide clamps	0.10/0.04 in.	103 in. (261cm)	0.3 mL (R) 0.3 mL (W) 0.3 mL (Y) 2.8 mL (B)*	50




*Priming volume does not include spike line.

To help clinicians accommodate an ever-changing patient population, Quest offers safe and simple solutions for fluid management. From high-flow stopcocks to the rigid MultiPort® Manifold and flexible ExoFlex™ Manifold, Quest designed the complete manifold line to adapt to your needs. The closed-system, procedure-ready sets aid in decreasing line manipulation and lowering infection risks, providing you even more time to devote to patient care.




95786	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (4) double T-sites with Q2® NCs + integrated check valves, removable spines and slide clamp	0.06 in.	11.25 in. (29cm)	1.75 mL	50




95787	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (6) T-sites with Q2® NCs + integrated check valves, removable spines and pinch clamp	0.06 in.	14.75 in. (37.5cm)	1.65 mL	50



95770	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	9 in. (22cm)	1.2 mL	50



95789	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	53 in. (135cm)	6.6 mL	50



95781	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	10.5 in. (26.75cm)	1.75 mL	50



95782	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	10.5 in. (26.75cm)	1.17 mL	50



95791	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	16 in. (40.65cm)	3 mL	50



95788	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	55 in. (140cm)	7.3 mL	50




95779	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	53 in. (135cm)	6.7 mL	50




95778	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	9 in. (22cm)	1.2 mL	50




95712	TUBING ID	APPROX. LENGTH	APV	CASE QTY
MultiPort® Manifold with (6) Q2® NCs + integrated check valves, double T-site with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	2.2 mL	50



95711	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (3) double T-sites with Q2® NCs + integrated check valves, (2) slide clamps and additional T-site with Q2® NC	0.06 in.	8.5 in. (22cm)	1.5 mL	50




9520	TUBING ID	APPROX. LENGTH	APV	CASE QTY
MultiPort® Manifold with (6) Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	2.0 mL	50



95710	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (6) T-sites with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	14.5 in. (37cm)	1.0 mL	50




95704	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves and slide clamp	0.06 in.	7 in. (18cm)	1.1 mL	50



95702	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	0.9 mL	50




95707	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs	0.06 in.	4 in. (10cm)	0.7 mL	50




High-Flow Stopcocks


ST01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
High-Flow 4-Way Stopcock	N/A	N/A	0.6 mL	100




ST02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
2-Gang High-Flow 4-Way Stopcock	N/A	N/A	0.8 mL	50




ST03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
3-Gang High-Flow 4-Way Stopcock	N/A	N/A	1.5 mL	50




STV01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	.5 mL	100



STV02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® 2-Gang High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	0.9 mL	50



STV03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® 3-Gang High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	1.1 mL	50



High-Flow Stopcocks continued


HFSTV01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® High-Flow 4-Way Stopcock Extension Set with Q2® NC, with pinch clamp + Y-site with Q2® NC	0.12 in.	43.8 in. (111cm)	7.7 mL	50



HFSTV02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® 2-Gang High-Flow 4-Way Stopcock Extension Set with Q2® NC, with pinch clamp + Y-site with Q2® NC	0.10 in.	45.5 in. (115cm)	8.0 mL	50



HFSTV03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® 3-Gang High-Flow 4-Way Stopcock Extension Set with Q2® NC, with pinch clamp + Y-site with Q2® NC	0.10 in.	47.1 in. (119cm)	8.3 mL	50




PRECISION DELIVERY® EXTENSION SETS




Quest Medical's Precision Delivery® Infusion Therapy Line is intentionally engineered to provide safe and ultra-precise infusions for your youngest critical-care patients. Specialized componentry minimizes dilution and priming volume while delivering all solutions to neonatal and pediatric patients. Reduce the risk of CLABSIs using our Q2® Valve Technology. A variety of configurations with specialized componentry such as ultra-microbore tubing, color-coded slide clamps, light-protective, amber tubing and pediatric-specific filters, provide for ultimate control of infusions.




PD06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 1.2 micron filter	0.04 in.	7.8 in. (19cm)	0.4 mL	50




PDM06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 0.2 micron filter	0.02 in.	7.8 in. (19cm)	0.3 mL	50




PDMV06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and 0.2 micron filter	0.02 in.	8 in. (20cm)	0.3 mL	50



PDV06	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC and slide clamp	0.04 in.	7.75 in. (19cm)	0.4 mL	50



PDV06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and 1.2 micron filter	0.04 in.	8.5 in. (21cm)	0.7 mL	50





PRECISION DELIVERY® EXTENSION SETS

PD60

Precision Delivery® Extension Set with slide clamp


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	61 in. (155cm)	1.4 mL	100



PD60F

Precision Delivery® Extension Set with slide clamp and 1.2 micron filter

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	62 in. (158cm)	1.7 mL	100



PDM60

Precision Delivery® Extension Set with slide clamp and ultra-microbore tubing

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	61 in. (155cm)	0.5 mL	100



PDM60F

Precision Delivery® Extension Set with slide clamp, 0.2 micron filter and ultra-microbore tubing

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	62 in. (158cm)	0.6 mL	100



PDMV60

Precision Delivery® Extension Set with Q2® NC, slide clamp and ultra-microbore tubing


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	61 in. (155cm)	0.6 mL	100



PDMV60F

Precision Delivery® Extension Set with Q2® NC, slide clamp, 0.2 micron filter and ultra-microbore tubing

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	62 in. (158cm)	0.8 mL	100



PDV60

Precision Delivery® Extension Set with Q2® NC and slide clamp

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	61 in. (155cm)	1.6 mL	100




PRECISION DELIVERY® EXTENSION SETS



PDV60F

Precision Delivery® Extension Set with Q2® NC, slide clamp and 1.2 micron filter


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	62 in. (158cm)	1.6 mL	100



PDBV05

Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps and 2-in-1 male luer


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	6 in. (16.5cm)	0.33 mL (W) 0.33 mL (B)	50



PDB05F

Precision Delivery® Bifuse Extension Set with (2) slide clamps, 0.2 micron filter and 2-in-1 male luer


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.28 in. (18.5cm)	0.38 mL (W) 0.55 mL (Y)	50



PDB05DF

Precision Delivery® Bifuse Extension Set with (2) slide clamps, 0.2 micron filter, 1.2 micron filter and 2-in-1 male luer


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.28 in. (18.5cm)	0.55 mL (W) 0.55 mL (Y)	50



PDBV05F

Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps, 0.2 micron filter and 2-in-1 male luer


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.4 in. (19cm)	0.33 mL (W) 0.50 mL (Y)	50



PDBV05DF

Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps, 0.2 micron filter, 1.2 micron filter and 2-in-1 male luer

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.4 in. (19cm)	0.50 mL (W) 0.50 mL (Y)	50






PRECISION DELIVERY® EXTENSION SETS

PDT05

Precision Delivery® Trifuse Extension Set with (3) slide clamps and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.2 mL (W)
0.2 mL (B)
0.2 mL (Y)

CASE QTY

50

PDT05F

Precision Delivery® Trifuse Extension Set with (3) slide clamps, 0.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.2 mL (W)
0.2 mL (B)
0.4 mL (Y)

CASE QTY

50

PDT05DF

Precision Delivery® Trifuse Extension Set with (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.4 mL (W)
0.2 mL (B)
0.4 mL (Y)

CASE QTY

50

PDTV05

Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.25 mL (W)
0.25 mL (B)
0.25 mL (Y)

CASE QTY

50

PDTV05F

Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps, 0.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.25 mL (W)
0.25 mL (B)
0.50 mL (Y)

CASE QTY

50

PDTV05DF

Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV

0.50 mL (W)
0.25 mL (B)
0.50 mL (Y)

CASE QTY

50

PRECISION DELIVERY® EXTENSION SETS




Light-Protective, Amber Extension Sets


PDV06A	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC and slide clamp	0.04 in.	7.75 in. (20cm)	0.4 mL	50




PDV60A	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC and slide clamp	0.04 in.	61 in. (155cm)	1.5 mL	50




PDMV60FA	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with ultra-microbore, light-protective, amber tubing, Q2® NC, slide clamp and 0.2 micron filter	0.02 in.	62 in. (158cm)	0.8 mL	50




PDV60FA	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC, slide clamp and 1.2 micron filter	0.04 in.	62 in. (158cm)	1.6 mL	50



PDTV05A	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.25 mL (Y)	50



PDTV05FA	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps, 0.2 micron filter and 3-in-1 male luer	0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.50 mL (Y)	50






PRECISION DELIVERY® EXTENSION SETS

Light-Protective, Amber Extension Sets continued

PDTV05DFA

Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV

0.50 mL (W)
0.25 mL (B)
0.50 mL (Y)

CASE QTY

50



EXTENSION SETS

Quest's Extension Sets focus on safe and accurate infusion therapy delivery with minimal priming volumes and high flow rates. Our offerings include various tubing sizes and component options including low-pressure, power injectable extensions that allow you to connect directly to your patient's catheter when administering contrast for diagnostic imaging.



95701	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Extension set with Q2® NC + integrated check valve, (1) T-site with Q2® NC + integrated check valve and pinch clamp	0.06 in.	7 in. (18cm)	0.66 mL	50




9219	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Female Y-extension set with (1) female luer pressure-activated check valve, (1) Q2® NC + integrated check valve and pinch clamp	0.10 in.	42 in. (107cm)	6.2 mL	50



95902	TUBING ID	APPROX. LENGTH	APV	CASE QTY
High-flow filter extension set with high-flow 0.22 micron filter, pinch clamp and Y-site with Q2® NC	0.10 in.	17 in. (43.2cm)	4.5 mL	100



95725S	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with 24 in. perforated sleeve and slide clamp	0.045 in.	25 ft. (7.6m)	10 mL	50



95724	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Extension set with pinch clamp	0.10 in.	20 ft. (6m)	31 mL	100



EXTENSION SETS


95725	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	20 ft. (6m)	8.5 mL	100



95726	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	10 ft. (3m)	3.2 mL	50



95727	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	60 in. (152cm)	1.7 mL	100




HF0135	TUBING ID	APPROX. LENGTH	APV	CASE QTY
High-flow extension set with (2) Y-sites + Q2® NCs and pinch clamp	0.100 in.	35 in. (88.9cm)	4.7 mL	50




Low-Pressure, Power Injection Extension Sets

Up to 325 psi


95714	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with detachable Q2®NC and pinch clamp	0.10 in.	7 in. (18cm)	1.0 mL	100



95906	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with slide clamp and Q2® NC + T-connector extension	0.045 in.	6 in. (15cm)	0.3 mL	100



95907	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with Q2® NC side port, slide clamp and Q2® NC + T-connector extension	0.045 in.	6 in. (15cm)	0.4 mL	100




Q2® NEEDLELESS CONNECTORS




Understanding patient safety is your #1 priority, Quest has revolutionized needleless connectors with our Q2® Valve Technology. Integrated into most Quest products, the enhanced safety features of Q2® NCs mitigate infection risks. The straight, unobstructed fluid path of Q2® NCs provides high flow rates and is engineered to suppress thrombotic occlusion risks caused by disconnecting syringes and other sets. Q2® NCs enable multiple activations and maintain a sterile, closed system for increased patient safety and improved outcomes, providing you with peace of mind.




3150	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® Needleless Connector, swabable	N/A	N/A	0.09 mL	200



3150HF	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Coming Soon Q2® High-Flow Needleless Connector, swabable	N/A	N/A	0.09 mL	200



3160	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® Vial Access Adapter, swabable	N/A	N/A	0.06 mL	50



QUEST Medical, Inc.

Ordering Process and Terms

To place an order for Infusion Therapy Products, please see below for ordering process and terms.

Distribution/Order Placement

Quest Medical, Inc.
One Allentown Parkway
Allen, Texas 75002

Phone: 800.627.0226 x 461

Fax: 972.390.7173

Email: Custserv@questmedical.com

Vendor Federal Tax I.D. No.

75-274056

F.O.B. ORIGIN

F.O.B., Pre-pay, and Add

Minimum Order

Need to order in case quantity

Delivery Lead Time

Lead time will be no more than forty-eight (48) to seventy-two (72) hours following receipt of order.

*Unless specified with 30-day lead time

Payment Terms

Net 30 days, with approval

Remit Address

Quest Medical, Inc.
PO Box 934486
Atlanta, GA 31193-4486

A/R Phone: 800.627.0226 x 271 or 299

***Custom Setup Inquiries** contact

Samples@questmedical.com

CONNECT WITH QUEST!





QUEST Medical, Inc.

custserv@questmedical.com | +1 800.627.0226
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com



EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Rx ONLY  Made in USA

© 2021 Quest Medical, Inc. Q2, ExoFlex, Precision Delivery, Qlinch and Multiport are registered trademarks of Quest Medical, Inc. 2021-20