




## Infusion Therapy Products



QUEST Medical, Inc.



**Quest Medical, Inc.** is committed to partnering with clinicians worldwide to produce better patient outcomes. We engineer infusion therapy products for superior infection control, organization, versatility and ease of use. Our innovative and high quality IV delivery sets are used to supply vital fluids and medications to patients in every department of your hospital or clinic.

Quest Medical aims to create efficiencies with each partnership through:

- Inventory Management**
- Reliable Client Services**
- Next Day Shipping**

Our representatives are subject matter experts who will assist you in standardizing products across departments, resulting in cost savings, reduced waste and faster set-up time.

To locate a sales specialist in your area, contact **800-672-0226** or email **[custserv@questmedical.com](mailto:custserv@questmedical.com)**.

#### Icon Description Key



Contains Q2\* Valve Technology



Low Volume (<3mL)



Contains Aspiration Port



European Conformity



Lead Time: 30 working days



#### BLOOD ADMINISTRATION SETS

Y115  
YV115  
YV130

5

#### IV ADMINISTRATION SETS

ASV10125  
ASV10115  
ASV10130  
ASV60125

6

ASV101Y  
ASV102Y  
ASV103Y  
95780  
22-201-V

7

#### MICROBORE ADMINISTRATION SETS

ASV170F  
ASV370  
ASV470

8

#### MANIFOLDS

95786  
95787  
95770  
95789

9

95781  
95782  
95791  
95788  
95779  
95778  
95712

10

95711  
9520  
95710  
95704  
95702  
95707

11

#### HIGH-FLOW STOPCOCKS

ST01  
ST02  
ST03  
STV01  
STV02  
STV03

12

HFSTV01  
HFSTV02  
HFSTV03

13

#### PRECISION DELIVERY® EXTENSION SETS

PD06F  
PDM06F  
PDMV06F  
PDV06  
PDV06F

14

PD60  
PD60F  
PDM60  
PDM60F  
PDMV60  
PDMV60F  
PDV60

15

PDV60F  
PDBV05  
PDB05F  
PDB05DF  
PDBV05F  
PDBV05DF

16

PDT05  
PDT05F  
PDT05DF  
PDTV05  
PDTV05F  
PDTV05DF

17

#### LIGHT-PROTECTIVE, AMBER EXTENSION SETS

PDV06A  
PDV60A  
PDMV60FA  
PDV60FA  
PDTV05A  
PDTV05FA

18

PDTV05DFA

19

#### EXTENSION SETS

95701  
9219  
95902  
95725S  
95724

20

95725  
95726  
95727  
HF0135

21

#### LOW-PRESSURE, POWER INJECTION EXTENSION SETS

95714  
95906  
95907

#### Q2® NEEDLELESS CONNECTORS

3150  
3150HF  
3160

22

## Our Evolving Product Lines


Quest Medical, Inc. was founded on the premise of developing medical devices that would have a positive impact on patient care. Since 1979, we have been dedicated to providing the highest quality products to our three focused markets: Infusion Therapy, Cardiovascular and Ophthalmology.

Quest's evolving Infusion Therapy Line features rigid and flexible manifolds with removable components and intravenous (IV) administration sets with high or low volume configurations to fit all patients' needs. Our products incorporate the following trademarks: Q2® Valve Technology (Needleless Connectors/NCs), Precision Delivery® Extension Sets, ExoFlex® Manifolds, MultiPort® Manifolds and Qlinch™ Securement Bands.



Integrated into most of Quest's Infusion Products is our Q2® Valve Technology. Gaining access to each patient's bloodstream via Infusion Therapy can increase risk of infection. Minimize that risk by incorporating a simple, closed-system design into your line setup. Quest's Q2® Valve Technology is void of internal componentry, reducing the chance of contamination.

**“When needleless systems are used, a split septum valve may be preferred over some mechanical valves due to increased risk of infection with the mechanical valves.”**




- CDC Guidelines for the Prevention of Intravascular Catheter-Related Infections; 2011:20, Department of Health & Human Services




## Blood Administration Sets


Y115	TUBING ID	APPROX. LENGTH	APV	CASE QTY
(2) high-flow vented spikes, (3) roller clamps, 200 micron filter chamber, large volume cylinder hand pump, slide clamp, (3) Y-sites + NCs, 32 in. detachable high flow 4-way 2-gang stopcock extension with pinch clamp, and detachable 6 in. extension with pinch clamp	.12 in.	119 in. (302 cm)	76 mL	50



YV115	TUBING ID	APPROX. LENGTH	APV	CASE QTY
(2) high-flow vented spikes, (3) roller clamps, 200 micron filter chamber, large volume cylinder hand pump, slide clamp and 35 in. detachable high flow 4-way 2-gang stopcock + Q2® NCs extension	.12 in.	114 in. (290 cm)	74 mL	50




YV130	TUBING ID	APPROX. LENGTH	APV	CASE QTY
(2) high-flow vented spikes, (3) roller clamps, 200 micron filter chamber, large volume cylinder hand pump, (2) slide clamps, (2) Y-sites + Q2® NCs, 45 in. detachable high flow 4-way 3-gang stopcock + Q2® NCs extension, and detachable 10 in. extension	.12 in.	129 in. (329 cm)	75 mL	50



# IV ADMINISTRATION SETS



Quest's IV Administration Sets assist with accurate delivery of medication. Accessing the patient's bloodstream when administering drugs may increase the risk of infection. With the incorporation of Q2® Needleless Connectors, Quest's IV Administration Sets are engineered with a closed system to promote patient safety. Multiple tubing ID options and detachable extensions provide a variety of configurations, which meet the demands of various departments in your facility including Anesthesia, Oncology and Critical Care.



## ASV10125


Vented 10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, 40 in. detachable high-flow 4-way 2-gang stopcock extension + Q2® NCs and pinch clamp

TUBING ID  
0.10 in.

APPROX. LENGTH  
129 in. (327cm)

APV  
28.6 mL

CASE QTY  
50



## ASV10115


Vented 10-drop/mL administration set with check valve, Y-site with Q2® NC, roller clamp, 22 in. detachable high-flow 4-way stopcock extension+ Q2® NCs and detachable 10 in. high-flow 4-way 2-gang stopcock extension + Q2® NCs

TUBING ID  
0.10 in.

APPROX. LENGTH  
115 in. (292cm)

APV  
27.9 mL

CASE QTY  
50



## ASV10130


Vented 10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, detachable 46 in. high-flow 4-way 3-gang stopcock extension + Q2® NCs and detachable 6 in. extension with Q2® NC and pinch clamp

TUBING ID  
0.10 in.

APPROX. LENGTH  
132 in. (335cm)

APV  
30.4 mL

CASE QTY  
50



## ASV60125


Vented 60-drop/mL administration set with check valve, (4) Y-sites with Q2® NCs, (2) pinch clamps, roller clamp, 35 in. detachable high-flow 4-way stopcock extension + Q2® NC and detachable 6.5 in. extension with Q2® NC

TUBING ID  
0.10 in.

APPROX. LENGTH  
126 in. (320cm)

APV  
29.1 mL

CASE QTY  
50






# IV ADMINISTRATION SETS


## ASV101Y

10-drop/mL administration set with check valve, (1) Y-site with Q2® NC, roller clamp, (1) T-site with Q2® NC + integrated check valve and pinch clamp




## ASV102Y

10-drop/mL administration set with check valve, (2) Y-sites with Q2® NCs, roller clamp, (1) T-site with Q2® NC + integrated check valve and pinch clamp




## ASV103Y

10-drop/mL administration set with check valve, (3) Y-sites with Q2® NCs, roller clamp, (2) T-sites with Q2® NCs + integrated check valve and pinch clamp




## 95780

10-drop/mL administration set with check valve, (3) Y-sites + Q2® NC, roller clamp, ExoFlex™ Manifold with (3) T-sites with Q2® NCs + integrated check valves + removable spine and detachable 45 in. extension set with pinch clamp



## 22-201-V

20-drop/mL vented secondary set with ball shut-off valve and roller clamp






## Microbore Administration Sets


### ASV170F

Microbore administration set with high-flow spike, pressure activated valve, pinch clamp, Q2® NC syringe port, check valve, slide clamp and 0.2 micron filter




### ASV370

Microbore trifurcated administration set with high-flow spike, pressure activated valve, pinch clamp, (3) Q2® NC syringe ports, check valve and (3) slide clamps



### ASV470

Microbore quadfurcated administration set with high-flow spike, pressure activated valve, pinch clamp, (4) Q2® NC syringe ports, check valve and (4) slide clamps




\*Priming volume does not include spike line.


To help clinicians accommodate an ever-changing patient population, Quest offers safe and simple solutions for fluid management. From high-flow stopcocks to the rigid MultiPort® Manifold and flexible ExoFlex™ Manifold, Quest designed the complete manifold line to adapt to your needs. The closed-system, procedure-ready sets aid in decreasing line manipulation and lowering infection risks, providing you even more time to devote to patient care.




95786	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (4) double T-sites with Q2® NCs + integrated check valves, removable spines and slide clamp	0.06 in.	11.25 in. (29cm)	1.75 mL	50




95787	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (6) T-sites with Q2® NCs + integrated check valves, removable spines and pinch clamp	0.06 in.	14.75 in. (37.5cm)	1.65 mL	50



95770	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	9 in. (22cm)	1.2 mL	50



95789	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	53 in. (135cm)	6.6 mL	50








95781	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	10.5 in. (26.75cm)	1.75 mL	50



95782	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	10.5 in. (26.75cm)	1.17 mL	50




95791	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	16 in. (40.65cm)	3 mL	50




95788	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	55 in. (140cm)	7.3 mL	50




95779	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs + integrated check valves, removable spine, pinch clamp and (1) Y-site with Q2® NC	0.10 in.	53 in. (135cm)	6.7 mL	50




95778	TUBING ID	APPROX. LENGTH	APV	CASE QTY
ExoFlex™ Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs + integrated check valves, removable spine and pinch clamp	0.10 in.	9 in. (22cm)	1.2 mL	50




95712	TUBING ID	APPROX. LENGTH	APV	CASE QTY
MultiPort® Manifold with (6) Q2® NCs + integrated check valves, double T-site with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	2.2 mL	50



95711	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (3) double T-sites with Q2® NCs + integrated check valves, (2) slide clamps and additional T-site with Q2® NC	0.06 in.	8.5 in. (22cm)	1.5 mL	50




9520	TUBING ID	APPROX. LENGTH	APV	CASE QTY
MultiPort® Manifold with (6) Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	2.0 mL	50



95710	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (6) T-sites with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	14.5 in. (37cm)	1.0 mL	50




95704	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (2) double T-sites with Q2® NCs + integrated check valves and slide clamp	0.06 in.	7 in. (18cm)	1.1 mL	50



95702	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (3) T-sites with Q2® NCs + integrated check valves and pinch clamp	0.06 in.	10 in. (25cm)	0.9 mL	50




95707	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Manifold with Q2® NC + integrated check valve, (2) T-sites with Q2® NCs	0.06 in.	4 in. (10cm)	0.7 mL	50




## High-Flow Stopcocks


ST01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
High-Flow 4-Way Stopcock	N/A	N/A	0.6 mL	100




ST02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
2-Gang High-Flow 4-Way Stopcock	N/A	N/A	0.8 mL	50




ST03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
3-Gang High-Flow 4-Way Stopcock	N/A	N/A	1.5 mL	50




STV01	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	.5 mL	100



STV02	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® 2-Gang High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	0.9 mL	50



STV03	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2® 3-Gang High-Flow 4-Way Stopcock with Q2® NC	N/A	N/A	1.1 mL	50



## High-Flow Stopcocks continued

### HFSTV01

Q2® High-Flow 4-Way Stopcock Extension Set with Q2® NC, with pinch clamp + Y-site with Q2® NC

TUBING ID

0.12 in.

APPROX. LENGTH

43.8 in. (111cm)

APV

7.7 mL

CASE QTY

50



### HFSTV02

Q2® 2-Gang High-Flow 4-Way Stopcock Extension Set with Q2® NC, with pinch clamp + Y-site with Q2® NC

TUBING ID

0.10 in.

APPROX. LENGTH

45.5 in. (115cm)

APV

8.0 mL

CASE QTY

50



### HFSTV03

Q2® 3-Gang High-Flow 4-Way Stopcock Extension Set with Q2® NC, with pinch clamp + Y-site with Q2® NC

TUBING ID

0.10 in.

APPROX. LENGTH



47.1 in. (119cm)

APV

8.3 mL

CASE QTY

50




# PRECISION DELIVERY® EXTENSION SETS




Quest Medical's Precision Delivery® Infusion Therapy Line is intentionally engineered to provide safe and ultra-precise infusions for your youngest critical-care patients. Specialized componentry minimizes dilution and priming volume while delivering all solutions to neonatal and pediatric patients. Reduce the risk of CLABSIs using our Q2® Valve Technology. A variety of configurations with specialized componentry such as ultra-microbore tubing, color-coded slide clamps, light-protective, amber tubing and pediatric-specific filters, provide for ultimate control of infusions.




PD06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 1.2 micron filter	0.04 in.	7.8 in. (19cm)	0.4 mL	50




PDM06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 0.2 micron filter	0.02 in.	7.8 in. (19cm)	0.3 mL	50




PDMV06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and 0.2 micron filter	0.02 in.	8 in. (20cm)	0.3 mL	50



PDV06	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC and slide clamp	0.04 in.	7.75 in. (19cm)	0.4 mL	50



PDV06F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and 1.2 micron filter	0.04 in.	8.5 in. (21cm)	0.7 mL	50








# PRECISION DELIVERY® EXTENSION SETS


PD60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp	0.04 in.	61 in. (155cm)	1.4 mL	100




PD60F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and 1.2 micron filter	0.04 in.	62 in. (158cm)	1.7 mL	100



PDM60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp and ultra-microbore tubing	0.02 in.	61 in. (155cm)	0.5 mL	100




PDM60F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with slide clamp, 0.2 micron filter and ultra-microbore tubing	0.02 in.	62 in. (158cm)	0.6 mL	100



PDMV60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp and ultra-microbore tubing	0.02 in.	61 in. (155cm)	0.6 mL	100



PDMV60F	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC, slide clamp, 0.2 micron filter and ultra-microbore tubing	0.02 in.	62 in. (158cm)	0.8 mL	100



PDV60	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Precision Delivery® Extension Set with Q2® NC and slide clamp	0.04 in.	61 in. (155cm)	1.6 mL	100




# PRECISION DELIVERY® EXTENSION SETS




## PDV60F

Precision Delivery® Extension Set with Q2® NC, slide clamp and 1.2 micron filter




## PDBV05

Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps and 2-in-1 male luer




## PDB05F

Precision Delivery® Bifuse Extension Set with (2) slide clamps, 0.2 micron filter and 2-in-1 male luer




## PDB05DF

Precision Delivery® Bifuse Extension Set with (2) slide clamps, 0.2 micron filter, 1.2 micron filter and 2-in-1 male luer




## PDBV05F

Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps, 0.2 micron filter and 2-in-1 male luer



## PDBV05DF

Precision Delivery® Bifuse Extension Set with (2) Q2® NCs, (2) slide clamps, 0.2 micron filter, 1.2 micron filter and 2-in-1 male luer






# PRECISION DELIVERY® EXTENSION SETS

## PDT05

Precision Delivery® Trifuse Extension Set with (3) slide clamps and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.2 mL (W)  
0.2 mL (B)  
0.2 mL (Y)

CASE QTY

50

## PDT05F

Precision Delivery® Trifuse Extension Set with (3) slide clamps, 0.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.2 mL (W)  
0.2 mL (B)  
0.4 mL (Y)

CASE QTY

50

## PDT05DF

Precision Delivery® Trifuse Extension Set with (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.4 mL (W)  
0.2 mL (B)  
0.4 mL (Y)

CASE QTY

50

## PDTV05

Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.25 mL (W)  
0.25 mL (B)  
0.25 mL (Y)

CASE QTY

50

## PDTV05F

Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps, 0.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV


0.25 mL (W)  
0.25 mL (B)  
0.50 mL (Y)

CASE QTY

50

## PDTV05DF

Precision Delivery® Trifuse Extension Set with (3) Q2® NCs, (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV

0.50 mL (W)  
0.25 mL (B)  
0.50 mL (Y)

CASE QTY

50



## Light-Protective, Amber Extension Sets

### PDV06A

Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC and slide clamp



### PDV60A


Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC and slide clamp



### PDMV60FA


*Coming Soon*

Precision Delivery® Light-Protective Amber Extension Set with ultra-microbore, light-protective, amber tubing, Q2® NC, slide clamp and 0.2 micron filter




### PDV60FA

Precision Delivery® Light-Protective Amber Extension Set with light-protective, amber tubing, Q2® NC, slide clamp and 1.2 micron filter




### PDTV05A

Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps and 3-in-1 male luer



### PDTV05FA

Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps, 0.2 micron filter and 3-in-1 male luer



TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.75 in. (20cm)	0.4 mL	50

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	61 in. (155cm)	1.5 mL	50

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.02 in.	62 in. (158cm)	0.8 mL	50

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	62 in. (158cm)	1.6 mL	50


TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.25 mL (Y)	50

TUBING ID	APPROX. LENGTH	APV	CASE QTY
0.04 in.	7.25 in. (18cm)	0.25 mL (W) 0.25 mL (B) 0.50 mL (Y)	50

## Light-Protective, Amber Extension Sets continued

### PDTV05DFA

Precision Delivery® Light-Protective Amber Trifuse Extension Set with light-protective, amber tubing, (3) Q2® NCs, (3) slide clamps, 0.2 micron filter, 1.2 micron filter and 3-in-1 male luer



TUBING ID

0.04 in.

APPROX. LENGTH

7.25 in. (18cm)

APV

0.50 mL (W)  
0.25 mL (B)  
0.50 mL (Y)

CASE QTY

50







Quest's Extension Sets focus on safe and accurate infusion therapy delivery with minimal priming volumes and high flow rates. Our offerings include various tubing sizes and component options including low-pressure, power injectable extensions that allow you to connect directly to your patient's catheter when administering contrast for diagnostic imaging.




95701	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Extension set with Q2® NC + integrated check valve, (1) T-site with Q2® NC + integrated check valve and pinch clamp	0.06 in.	7 in. (18cm)	0.66 mL	50




9219	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Female Y-extension set with (1) female luer pressure-activated check valve, (1) Q2® NC + integrated check valve and pinch clamp	0.10 in.	42 in. (107cm)	6.2 mL	50



95902	TUBING ID	APPROX. LENGTH	APV	CASE QTY
High-flow filter extension set with high-flow 0.22 micron filter, pinch clamp and Y-site with Q2® NC	0.10 in.	17 in. (43.2cm)	4.5 mL	100



95725S	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with 24 in. perforated sleeve and slide clamp	0.045 in.	25 ft. (7.6m)	10 mL	50



95724	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Extension set with pinch clamp	0.10 in.	20 ft. (6m)	31 mL	100




95725	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	20 ft. (6m)	8.5 mL	100



95726	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	10 ft. (3m)	3.2 mL	50



95727	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Microbore extension set with slide clamp	0.045 in.	60 in. (152cm)	1.7 mL	100




HF0135	TUBING ID	APPROX. LENGTH	APV	CASE QTY
High-flow extension set with (2) Y-sites + Q2® NCs and pinch clamp	0.100 in.	35 in. (88.9cm)	4.7 mL	50




## Low-Pressure, Power Injection Extension Sets

Up to 325 psi


95714	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with detachable Q2®NC and pinch clamp	0.10 in.	7 in. (18cm)	1.0 mL	100



95906	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with slide clamp and Q2® NC + T-connector extension	0.045 in.	6 in. (15cm)	0.3 mL	100



95907	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Low-pressure, power injection extension set with Q2® NC side port, slide clamp and Q2® NC + T-connector extension	0.045 in.	6 in. (15cm)	0.4 mL	100







Understanding patient safety is your #1 priority, Quest has revolutionized needleless connectors with our Q2<sup>®</sup> Valve Technology. Integrated into most Quest products, the enhanced safety features of Q2<sup>®</sup> NCs mitigate infection risks. The straight, unobstructed fluid path of Q2<sup>®</sup> NCs provides high flow rates and is engineered to suppress thrombotic occlusion risks caused by disconnecting syringes and other sets. Q2<sup>®</sup> NCs enable multiple activations and maintain a sterile, closed system for increased patient safety and improved outcomes, providing you with peace of mind.




3150	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2 <sup>®</sup> Needleless Connector, swabable	N/A	N/A	0.09 mL	200



3150HF	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2 <sup>®</sup> High-Flow Needleless Connector, swabable	N/A	N/A	0.09 mL	200



3160	TUBING ID	APPROX. LENGTH	APV	CASE QTY
Q2 <sup>®</sup> Vial Access Adapter, swabable	N/A	N/A	0.06 mL	50



# QUEST Medical, Inc.

## Ordering Process and Terms

To place an order for Infusion Therapy Products, please see below for ordering process and terms.

### Distribution/Order Placement

Quest Medical, Inc.  
One Allentown Parkway  
Allen, Texas 75002

**Phone:** 800.627.0226 x 461

**Fax:** 972.390.7173

**Email:** Custserv@questmedical.com

### Vendor Federal Tax I.D. No.

75-274056

### F.O.B. ORIGIN

F.O.B., Pre-pay, and Add

### Minimum Order

Need to order in case quantity

### Delivery Lead Time

Lead time will be no more than forty-eight (48) to seventy-two (72) hours following receipt of order.

\*Unless specified with 30-day lead time

### Payment Terms

Net 30 days, with approval

### Remit Address

Quest Medical, Inc.  
PO Box 934486  
Atlanta, GA 31193-4486

**A/R Phone:** 800.627.0226 x 271 or 299

**\*Custom Setup Inquiries** contact

Samples@questmedical.com

## CONNECT WITH QUEST!





**QUEST** Medical, Inc.

custserv@questmedical.com | +1 800.627.0226  
+1 972.390.9800 | fax: +1 972.390.7173 | questmedical.com



**EC REP**

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

**Rx ONLY**  Made in USA